

INSTRUMENT PROTECTION TRAYS

Your Trusted Leader for Instrument Longevity.



Innovia Medical is a **family of specialty surgical companies** that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

ELEVATING CARE. IMPROVING OUTCOMES.

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

INNOVATION LEADERSHIP.

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

UNSURPASSED PRODUCT QUALITY.

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.


From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.

We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at **www.innoviamedical.com**.

THE INSTRUSAFE ADVANTAGE	
360° OF PROTECTION	4
INDUSTRY DISADVANTAGES	6
FEATURES + BENEFITS	8
ORDER PROCESS: THE INSTRUSAFE WAY	10
STERILIZATION BARRIERS + FOOT OPTIONS	12
TRAYS BY PROCEDURE	
ARTHROSCOPY	14
BARIATRIC	16
CARDIOVASCULAR	18
DENTAL	20
ENT	22
GENERAL	34
GYNECOLOGY	38
NEUROLOGY/SPINE	40
OPHTHALMOLOGY	46
ORTHOPEDIC	50
PLASTICS	52
RETRACTORS	54
ROBOTICS	58
SCOPES	79
	77777

TRANSPORT CONTAINERS	
XL TRANSPORT CONTAINERS	90
HEAVY DUTY TRANSPORT CONTAINERS	92
TRANSPORT CONTAINERS	93
	////////
ACCESSORIES	
SILICONE SMALL PARTS BOXES	94
PLAQUE OPTIONS	95
FINGER MATS	96
GUARDIAN TRAY	96
INSTRU-ZYME	97
ROUNDERS POINT GUARDS	97
	////////



SUBMIT PURCHASE ORDERS TO:

Summit Medical LLC 815 Vikings Parkway, Suite 100 St. Paul, MN 55121 | USA

P: 1-888-229-2875 | +1 651-789-3939 F: 1-888-229-1941 | +1 651-789-3979 www.summitmedicalusa.com



The InstruSafe® Advantage

FDA 510(k) Clearance

Designed to ensure optimal sterilization, InstruSafe Trays have undergone rigorous testing before receiving multiple FDA 510(k) clearances. Cleared for various sterilization cycles, these reliable trays give hospitals and healthcare facilities even more options when it comes to organizing, sterilizing, and protecting their instruments.

For the most up-to-date list of cleared cycles, visit www.instrusafe.com/validations.




POINT CAMERA TO SEE VALIDATIONS

360° of Protection

The use of complex and delicate instrument sets is a new reality in which sterile processing and the operating room live. Intricacy of surgical instrumentation continues to increase, and so does the cleaning and sterilization process. Now, healthcare facilities face the challenge of reprocessing instruments safely and efficiently while also meeting the demand for fast turnaround. Keeping these sets organized and protected upon arrival in the operating room creates a challenge for all involved. Relieve yourself of these daily obstacles with InstruSafe Trays—an everyday solution for the intricacies of your instrument sets.

We have spent over 30 years understanding the use cycle and improving outcomes. With the use cycle in mind, we have created a product that eliminates common frustrations associated with instrument sets. The products in this catalog are a compilation of our hard work. As you journey through the pages, let us help you solve your everyday instrument problems

Summit Medical's most recent FDA 510(k) clearances and product advancements are a testament of our dedication to our customers' satisfaction. We wouldn't have continuing success without the healthcare community we serve. Focus, determination and your support will continue our impact on patient care for decades to come.



DECONTAMINATION

Instruments are individually protected going into the wash cycle and the trays are designed with large perforations, allowing maximum detergent and water flow.

TRAY PREPARATION

Each tray is designed specifically for the instruments within it, providing individual locations for each one. This protects them from harming the technician, from damaging themselves and from puncturing or jeopardizing the integrity of the sterile barrier.

STERILIZATION

The large amount of perforations in our trays provide superior sterilant penetration. Many of our silicone parts are also designed with minimal contact points and additional ventilation in each slot that match our tray holes to improve sterilant penetration to each individual instrument.

TRANSPORTATION

Soft and durable silicone holders in the base and cover secure instruments to prevent shifting and provide 360 degrees of protection. InstruSafe Trays help to avoid unnecessary damage and to present an organized instrument set every time.

PERATING ROOM USE

Each instrument is easily identifiable and placed in a singular slot for surgical presentation and tracking at the end of a procedure.



Designed to organize and protect from the point of operating room use to decontamination, our solutions benefit everyone who handles our products. THE INSTRUSAFE ADVANTAGE

THE INSTRUSAFE ADVANTAGE

Look Familiar?

Unfortunately, this is what most sterile processing departments have come to know. The familiarity and acceptance of these messy, hazardous baskets makes us cringe! Why, you ask? Because the minimal organization and protection of surgical instruments practiced in many facilities presents numerous risks to the individuals that come in contact with them, as well as to the instruments themselves.

DISADVANTAGES OF MESSY BASKETS

POTENTIAL TO COMPROMISE STERILE BARRIER BY POKING THROUGH WRAP -

DANGEROUS TO TRANSPORT FROM THE O.R. TO THE SPD -

MINIMAL INSTRUMENT PROTECTION AND ORGANIZATION -

DIFFICULT TO COUNT AND TRACK INSTRUMENTS -

BASKETS NEED FREQUENT REPLACEMENT

ADDS TIME TO SURGICAL SET-UP -

BROKEN INSTRUMENTS IF DROPPED -

HIGH RISK FOR FINGER POKES □

POOR SURGICAL PRESENTATION -

HARD TO FIND WHAT YOU NEED -

MISSING/LOST INSTRUMENTS -

UNNECESSARY REPROCESSING -


DAMAGED INSTRUMENTS -

WASTE OF SPACE -

EXPECT MORE.


For these reasons and many more, InstruSafe Trays offer custom solutions to meet each customer's specific needs. We think it's time for everyone to expect more from their trays. THE INSTRUSAFE ADVANTAGE

THE INSTRUSAFE ADVANTAGE



Simple Steps

If one of our preconfigured tray designs does not fit your needs let our product experts design a custom tray with our effortless process below.





ASSESS

An InstruSafe
Representative
conducts a free,
on-site assessment
of your instrument
set.



GATHER

Your instrument set specifics and inputs are gathered by our team of experts.



DESIGN

Using your inputs, our product specialists build a tray design with you and your instrument set in mind.



EVALUATE

You be the judge! Evaluate the new tray design while it's loaded with your instrument set. Let us know your thoughts and we will make any and all adjustments necessary!



COMPLETE

It's all yours! Once complete, each custom tray is assigned its own unique part number, making it simple and easy for you to reorder in future.

HOW TO USE THE CATALOG


The colored headers identify which procedure the instrument set you're looking at is used for. The trays on each page are differentiated by preconfigured or custom configured with the symbols shown in the key below. Some trays have several layers including a cover, shelves, baskets, and a base. The custom configured tray examples included in the catalog are a small representation of the trays our local representatives have built for facilities around the world.

CUSTOM + PRECONFIGURED SYMBOLS



THIS SYMBOL IDENTIFIES A PRECONFIGURED TRAY

Order by part number from your local InstruSafe representative or call us today.



THIS SYMBOL IDENTIFIES A

CUSTOM CONFIGURED TRAY EXAMPLE
Call us today to schedule your free
assessment.



Sterilization Barriers


InstruSafe Instrument Protection Trays are FDA cleared for several sterilization cycles using wrap and/or rigid containers.

For the most up-to-date list of FDA 510(k) cleared cycles, visit www.instrusafe.com.

Foot Options

R Foot: Rounders Corner Guards—to be used with wrap

TF Foot: Tall Foot fits into rigid containers



HOW TO ORDER:

EXAMPLE: TRAY IN-2890

Wrapped Foot Container Foot IN-2890-R IN-2890-TF

FOOT OPTION IS ADDED TO THE END OF THE PART NUMBER

Rounders® Corner Guards

Rounded, medical grade corner guards affix to InstruSafe Tray corners to reduce potential damage to the sterile wrap barrier. Rounders are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to all InstruSafe Trays.



ROUNDERS CORNER GUARDS

Small Rounders	
Kit of Four	IN-7334
Large Rounders	
Kit of Four	IN-7337

- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit all InstruSafe Trays
- · Rounded, medical grade silicone corner guard
- Help avoid compromised sterility in the O.R.
- Environmentally and economically efficient
- Impact resistant
- US Patent #9119611 | #D783856

Tall Foot

Designed to affix to trays that are intended to be sterilized within a rigid container. Tall Feet are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to most InstruSafe Trays.



TF FOOT KIT

Kit of Four	IN-7338

- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit most InstruSafe Trays

ARTHROSCOPY ARTHROSCOPY




over

Shelf

Preconfigured Tray for

Arthroscopy

This three-level tray is equipped for both small and large arthroscopic instrument sets. Our tray is sized to fit most threequarter size containers or can be wrapped for your convenience.



Custom Configured Tray Example for

Arthroscopy Example

ARTHROSCOPY INSTRUMENT TRAY

Arthroscopy Instrument Tray

A combination of up to 12 of the following instruments:

Rongeurs

Graspers

- Ring Handle Instruments
- Other Associated Endoscopic
- Instruments


15 ¼ L x 9 5% W x 5 3¼ H (in)	387 L x 244 W x 146 H (mm)
Wrapped Foot	Container Foot
IN-2846-R	IN-2846-TF

Expand capacity with positionable Hold-Its and Hold-Downs™ silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.

Basket and base can be customized to accommodate specific instrumentation. ARTHROSCOPY TMJ TRAY | 17 INSTRUMENTS






Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

BARIATRIC BARIATRIC



Preconfigured Tray for Bariatric

The bariatric tray is designed to accommodate even the longest instruments. Our three-level system allows you to get organized and protect all of your bariatric instrumentation in one tray.



Custom Configured Tray Example for

Bariatric

BARIATRIC TRAY

Bariatric Instrument Tray

A combination of up to 15 (32cm to 45cm+) of the following instruments:

- Laparoscopic Instruments • General Dissecting Instruments
- Cannulas
- Miscellaneous Parts

28 L x 13 W x 8 H (in)

- Sheaths • Trocars
- Tubing/Cords

Wrapped Foot IN-8864-R

Container Foot IN-8864-TF


711 L x 330 W x 203 H (mm)

Expand capacity with positionable Hold-Its and Hold-Downs silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.

CALL US TODAY FOR CUSTOMIZATION!

CUSTOM CONFIGURED BARIATRIC TRAY | 31 INSTRUMENTS



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

BARIATRIC LAPAROSCOPIC SET | 59 INSTRUMENTS


BEFORE





Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Cardiovascular




EVANS ROTATABLE TRAY | 32 INSTRUMENTS

Custom Configured Tray Example for


Cardiovascular

MICROVASCULAR SET | 19 INSTRUMENTS






CUSTOM CONFIGURED



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection



Dental


Custom Configured Tray Example for

Dental

DENTAL SET | 34 INSTRUMENTS

.... PEEODE






Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

OPERATIVE TRAY SET | 11 INSTRUMENTS

CUSTOM CONFIGURED



ARTHRO

GYN

NEURO

EYE

PLASTICS

SCORES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection

ENT

ENTSYSTEMS

Since the mid-1980's, we have offered a complete solution for fine and delicate ear instrumentation

through our Middle Ear System. We are proud to claim this as our very first InstruSafe product.

Our modules are designed to specifically fit each instrument and can be used individually or in

MIX AND MATCH OUR MODULES FROM PAGES 23-24 TO FIT YOUR

combination, to organize and protect an entire instrument set.

INSTRUMENT SET, CREATING YOUR UNIQUE ENT SYSTEM!



Preconfigured Trays for

ENT Modules







19 5% L x 9 5% W x 6 H (in)


Wrapped Foot

IN-2890-R

EAR MODULE TRAY

• 1 Speculum Module • 9½ x 13" Finger Mat 498 L x 244 W x 152 H (mm) Container Foot

IN-2890-TF





Alligator Module • 8 Alligator Style Instruments (Cup, Scissors, etc.) 5 % L x 5 % W x 1 H (in) 149 L x 130 W x 25 H (mm) IN-2001



EAR MODULE TRAY Pick + Hook Module, Wide • 10 Small Handle Picks, Hooks, Dissectors or Curettes 7 % L x 5 % W x 1 H (in) 194 L x 130 W x 25 H (mm) IN-2002



EAR MODULE TRAY

Speculum Module Holds: • 6-12 Speculas 3 ½ L x 5 ½ W x 1 H (in) 89 L x 130 W x 25 H (mm) IN-2003

1.888.229.2875 www.innoviamedical.com/instrusafe

Trav shown without cover for illustrative purposes

ENT



Preconfigured Trays for


ENT Modules



Preconfigured Trays for Picks + Hooks

Similar to our original module designs, these standalone trays will protect and organize pick, hook and knife instrumentation for ENT procedures. Instruments can either be positioned at a 45 degree angle for procedural presentation or secured flat for sterilization and storage.





EAR MODULE TRAY

Micro-Suction Module

• 12 Micro-Suction Irrigators with Luer Lock Hubs

5 ½ L x 3 ½ W x 1 H (in)

130 L x 89 W x 25 H (mm)

IN-2004



Suction Module

Holds:

• 10 Ferguson Style Suctions

5 % L x 5 % W x 1 H (in)

149 L x 130 W x 25 H (mm)

IN-2005



Alligator Module With Cover

• 8 Alligator Style Instruments (Cup, Scissors, etc.)

7 % L x 5 % W x 3 H (in)

194 L x 130 W x 76 H (mm)

IN-2006

EAR MODULE TRAY

Pick + Hook Module

Holds:

• 9 Small Handle Picks, Hooks, Dissectors or Curettes

194 L x 117 W x 25 H (mm) 7 % L x 4 % W x 1 H (in)

IN-2008

PICK + HOOK TRAY

10 Pick + Hook Tray

• 10 Pick + Hook Style Instruments

5 ½ L x 8 ¾ W x 1 ½ H (in) 133 L x 222 W x 38 H (mm)

IN-4010



PICK + HOOK TRAY

20 Pick + Hook Tray

• 20 Pick + Hook Style Instruments

11 L x 8 ³/₄ W x 1 ½ H (in)

279 L x 222 W x 38 H (mm)

IN-4000

PICK + HOOK TRAY

30 Pick + Hook Tray


• 30 Pick + Hook Style Instruments

15 L x 8 3/4 W x 1 1/2 H (in)

381 L x 222 W x 38 H (mm)

IN-4003

ENT ENT




Custom Configured Tray Example for


Ear

Custom Configured Tray Example for

Ear


MYRINGOTOMY SET | 26 INSTRUMENTS






Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.


EAR TRAY | 86 INSTRUMENTS





Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection

ENT ENT




Custom Configured Tray Example for


Ear



Preconfigured Tray for Nose/Sinus

MAJOR EAR TRAY | 66 INSTRUMENTS







Basket and base can be customized to accommodate specific instrumentation.

ENDOSCOPY SINUS TRAY

Endoscopy Sinus Instrument Tray

A combination of up to 12 of the following instruments:

Punches

 Rongeurs Graspers

• Ring Handle Instruments • Other Associated Endoscopic

IN-2846-R


Instruments 15 ½ L x 9 5/8 W x 5 3/4 H (in) 387 L x 244 W x 146 H (mm)

IN-2846-TF

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Summit Medical An Innovia Medical Company

NT ENT




Custom Configured Tray Example for

Nose/Sinus



SINUS SET | 32 INSTRUMENTS

BEFORE



SINUS TRAY | 34 INSTRUMENTS

CUSTOM CONFIGURED **Basket**

Base

TRANSPORT CONTAINER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

1.888.229.2875 www.innoviamedical.com/instrusafe

3

222 227

ENT ENT




Preconfigured Tray for Throat


Custom Configured Tray Example for

Throat

THYROID SET | 25 INSTRUMENTS



BEFORE



TRANSPORT CONTAINERS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protec


GENERAL GENERAL



Preconfigured Tray for


General

Maximize the protection and organization of your laparoscopic instrument sets with this three-level tray that is excellently constructed to fit both small and large sets of instrumentation.




Custom Configured Tray Example for

General



Basket and base can be customized to accommodate specific instrumentation. **TENDON STRIPPER | 13 INSTRUMENTS**

BEFORE




Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Base

1.888.229.2875 www.innoviamedical.com/instrusafe

General



Custom Configured Tray Example for

General


LAPAROSCOPIC | 47 INSTRUMENTS


*CANTON OF THE PARTY OF THE PAR CUSTOM -----CONFIGURED

Basket


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

BASIC UROLOGY TRAY | 35 INSTRUMENTS






Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Hysteroscopy




Custom Configured Tray Example for

Hysteroscopy

HYSTEROSCOPE SET | 10 INSTRUMENTS

Tray shown without cover for illustrative purposes, All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

BEFORE



HYSTERECTOMY RESECTOSCOPE TRAY | 15 INSTRUMENTS

CUSTOM CONFIGURED



Tray shown without cover for illustrative purposes, All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

NEUROLOGY/SPINE NEUROLOGY/SPINE



Preconfigured Trays for **Kerrisons**

Our Neurology and Spine Trays provide maximum accessibility and flexibility when handling and sterilizing kerrisons and forceps.



Preconfigured Trays for

Bi-Polar Forceps



BARIATRIC

CARDIC

DENTAL

Tray shown without cover for illustrative purposes.

KERRISON INSTRUMENT TRAY + RACK

Ten Position Instrument Tray

Holds:

• 10 Kerrisons

15 ½ L x 9 ½ W x 6 H (in)

Wrapped Foot
IN-5012-R

Container Foot
IN-5012-TF



BAYONET FORCEPS MODULE TRAY

Dual Forceps Module Tray with Cover Holds:

• 2 Forceps

9 ³/₄ L x 2 ½ W x 1 ½ H (in) 248 L x 54 W x 29 H (mm)

IN-5401-02

ENT

GENER

GYI

NEU

ROBOTIC

TRANSPORT CONTAINERS



KERRISON INSTRUMENT RACK

KERRISON INSTRUMENT RACK

Five Position Instrument Rack


• 5 Kerrisons

14 ½ L x 5 W x 6 H (in)

Wrapped Foot

N/A

Ten Position Instrument Rack	
Holds:	
• 10 Kerrisons	
14 ¼ L x 9 ½ W x 6 H (in)	362 L x 241 W x 152 H (mm)
Wrapped Foot N/A	Container Foot IN-5013-TF



BAYONET FORCEPS MODULE TRAY

Dual Extra Long Forceps Module Tray with Cover Holds:

• 2 Forceps

12 L x 2 1/8 W x 1 1/2 H (in) 305 L x 54 W x 38 H (mm)

IN-5401-03

101 In the second secon

BAYONET FORCEPS MODULE TRAY

Six Forceps Module Tray with Cover Holds:

• 6 Forceps

9 ³/₄ L x 4 ⁵/₈ W x 1 ¹/₈ H (in)

248 L x 248 W x 29 H (mm)

IN-5401

TRANSPORT CONTAINER

ACCESSORIES

ACCESSORIE

362 L x 127 W x 152 H (mm)


IN-5014-TF

NEUROLOGY/SPINE NEUROLOGY/SPINE



Preconfigured Trays for

Bi-Polar Forceps



Custom Configured Tray Example for

Neurology/Spine



BAYONET FORCEPS MODULE TRAY

Eight Forceps Module Tray with Cover

• 8 Forceps

305 L x 184 W x 38 H (mm) 12 L x 7 ¼ W x 1 ½ H (in)

IN-5401-08



BAYONET FORCEPS MODULE TRAY

Twelve Forceps Module Tray with Cover

• 12 Forceps

12 L x 9 5/8 W x 1 1/2 H (in)

305 L x 244 W x 38 H (mm)

IN-5401-12


NEUROLOGY TRAY | 13 INSTRUMENTS

BEFORE






Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




Neurology/Spine




Custom Configured Tray Example for

Neurology/Spine

HIGH SPEED PNEUMATIC DRILL TRAY | 9 INSTRUMENTS



DURAL REPAIR KIT TRAY | 18 INSTRUMENTS




Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.


1.888.229.2875

Summit Medical An Innovia Medical Company

InstruSafe Sterile Processing




Ophthalmology




Custom Configured Tray Example for

Ophthalmology

CATARACT-PHACO SET | 28 INSTRUMENTS








EYE TRAY | 57 INSTRUMENTS







Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection

1.888.229.2875 **Summit Medical** InstruSafe 46 www.innoviamedical.com/instrusafe Sterile Processing



Ophthalmology




Custom Configured Tray Example for

Ophthalmology


CORNEA TRAY | 61 INSTRUMENTS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

CUSTOM CONFIGURED



PHACO HAND PIECE TRAY | 23 INSTRUMENTS



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection

InstruSafe


Sterile Processing

ORTHOPEDIC ORTHOPEDIC



Preconfigured Trays for Orthopedic

While not usually delicate, orthopedic instrumentation needs to be kept sharp and straight to get the job done. Our Ortho Trays are built to safely accommodate and protect instruments that are commonly found in almost all surgical facilities.



Custom Configured Tray Example for

Orthopedic

ORTHOPEDIC TRAY

Lambotte Osteotome Tray

• 12 Instruments—straight or curved, up to 1 $\frac{1}{2}$ " width (38mm)

8 3% L x 9 5% W x 3 H (in)	213 L x 244 W x 76 H (mm)
Wrapped Foot	Container Foot
IN-5010-R	IN-5010-TF


CURETTE TRAY


Curette Tray

• 14 Instruments-7-13" in length


15 ¼ L x 9 5/8 W x 2 H (in)	387 L x 244 W x 51 H (mm)
Wrapped Foot	Container Foot
IN-5011-R	IN-5011-TF

SMALL FRAG SET | 23 INSTRUMENTS






Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Plastics




Custom Configured Tray Example for

Plastics

MICRO SET | 32 INSTRUMENTS

BEFORE



PLASTIC TRAY | 38 INSTRUMENTS

CUSTOM **CONFIGURED**

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

RETRACTORS RETRACTORS



Preconfigured Trays for Bookwalter® Retractors

This tray system is engineered to have a balanced weight distribution that meets the AAMI 25 lb. weight limit standard. We offer a two-level tray design made to accommodate the ring and blade components with an additional tray to store the table post and horizontal bar. Our two-tray system is organized for easy counting and safe for transporting. *Note:* The Mediflex Bookler Retractor also fits in this tray.



Custom Configured Tray Example for

Retractors

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex

BOOKWALTER RETRACTOR TRAY

Ring + Blade Tray

Holds:

- 1 (1") Post Coupling
- 1 Balfour Retractor • 1 Horizontal Flex Bar
- 1 Harrington Retractor
- 1 Medium Oval Ring
- 1 Segmented Ring
- 1 Gelpi Retractor
- 2 (2x6"/51mm x 152mm) Malleable Retractors
- 1 (3x6"/76mm x 152mm) Malleable Retractor
- 2 (2x3"/51mm x 76mm) Kelly Retractors
- 1 (2x4"/51mm x 110mm) Kelly Retractor
- 1 (2x5"/51mm x 127mm) Kelly Retractor
- 1 (2x6"/51mm x 152mm) Kelly Retractor
- 6 Tilt Ratchet Mechanisms
- 2 Ratchet Mechanisms

22 ½ L x 13 W x 4 ½ H (in)	572 L x 330 W x 114 H (mm)
Wrapped Foot	Container Foot
IN-8700-R	IN-8700-TF

BOOKWALTER RETRACTOR TRAY


Post Tray


Holds:

• Table Post and Horizontal Bar

	. 5
25 L x 5 3/8 W x 3 H (in)	635 L x 137 W x 76 H (mm)
Wrapped Foot	Container Foot
IN-8701-R	IN-8701-TF

GREENBERG SET | 17 INSTRUMENTS






Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

1.888.229.2875


Summit Medical An Innovia Medical Company InstruSafe Sterile Processing

1.888.229.2875 www.innoviamedical.com/instrusafe

www.innoviamedical.com/instrusafe




Retractors



Custom Configured Tray Example for

Retractors


BOOKWALTER RETRACTOR TRAY | 38 INSTRUMENTS



Shelf



CARDIAC MIS RETRACTOR HOLDER TRAY | 2 INSTRUMENTS



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing 1.888.229.2875 57




Preconfigured Trays for da Vinci SP®

The storage and reprocessing of the da Vinci SP isn't limited to a single solution. That's why we've designed multiple trays specific to the system's instrumentation to accommodate a wide variety of customer preferences.



Preconfigured Trays for

da Vinci SP




INSTRUMENT TRAY

Single Instrument Tray Holds: • 1 Instrument 679 L x 175 W x 95 H (mm) 26 ³/₄ L x 6 ⁷/₈ W x 3 ³/₄ H (in) Wrapped Foot Container Foot IN-8921-R IN-8921-TF



INSTRUMENT TRAY

Six Instrument Tray Holds: • 6 Instruments 711 L x 279 W x 137 H (mm) 28 L x 11 W x 5 3/8 H (in) Wrapped Foot Container Foot IN-8925-R IN-8925-TF



INSTRUMENT TRAY

INSTRUMENT TRAY

Four Instrument Tray

28 L x 8 W x 5 3/8 H (in)

Wrapped Foot

IN-8923-R

4 Instruments

Dual Instrument Tray	
Holds:	
• 2 Instruments	
26 ³ / ₄ L x 6 ⁷ / ₈ W x 3 ³ / ₄ H (in)	679 L x 175 W x 95 H (mm)
Wrapped Foot	Container Foot
IN-8922-R	IN-8922-TF



PROCEDURE TRAY

Four Instrument Procedure Tray Holds: • 4 Instruments + Associated Accessories 28 L x 8 W x 5 3/8 H (in) 711 L x 203 W x 137 H (mm) Wrapped Foot Container Foot IN-8924-R IN-8924-TF

PROCEDURE TRAY

Six Instrument Procedure Tray • 6 Instruments + Associated Accessories 28 L x 11 W x 5 3/8 H (in) 711 L x 279 W x 137 H (mm) Wrapped Foot Container Foot IN-8926-R IN-8926-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company

711 L x 203 W x 137 H (mm)

Container Foot

IN-8923-TF

InstruSafe Sterile Processing



Preconfigured Trays for

da Vinci SP



da Vinci Xi® Single-Site®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Single-Site instruments—from single instruments to procedural sets and their accommodating accessories—our trays are the recommended solution for instrument reprocessing.

ARTHRO

BARIATRI

CARDIO

DENTA

NEU

DI ASTI

SCOPES

TRANSPORT CONTAINERS



□ ACCESSORY TRAY

Six Accessory Tray

Holds:

2 Cannulas 2 Obturators

• 2 Cables

CAMERA TRAY

Holds:
• 1 Camera

Single Camera Tray


28 L x 11 W x 3 H (in)

Wrapped Foot

IN-8920-R

15 ¼ L x 9 5% W x 2 3% H (in) 387 L x 244 W x 60 H (mm)

Wrapped Foot Container Foot
IN-8927-R IN-8927-TF





□ INSTRUMENT TRAY

Single Instrument Tray

Holds:

• 1 (5mm) Xi Single-Site Instrument

27 ½ L x 3 ¾ W x 3 ½ H (in) 689 L x 86 W x 89 H (mm)

Wrapped Foot Container Foot N/A

PROCEDURE TRAY

22 Instrument Procedure Tray

olds:

• 6 (5mm) Xi Single-Site Instruments

• 3 Rigid Obturators

• 2 Flexible Obturators

• 4 Curved Cannulas

• 3 Rigid Cannulas

• 3 Cables

• 1 Suction Irrigator Tube and Connector

26 ¾ L x 13 W x 3 5% H (in) 679 L x 330 W x 92 H (mm)

Wrapped Foot
IN-8912-R

Container Foot
N/A

ACCESSORY TRAY

16 Accessory Tray

Holds:

3 Rigid Obturators 2 Flexible Obturators

• 3 Rigid Cannulas

• 4 Curved Cannulas

• 3 Cables

• 1 Suction Irrigator Tube and Connector

15 ¼ L x 9 % W x 4 % H (in)

Wrapped Foot
IN-8913-R

387 L x 244 W x 105 H (mm)

Container Foot
IN-8913-TF

All da Vinci Xi Single-Site Trays have green covers.

Excludes IN-8911

TRANSPORT CONTAINE

050000150

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ACCESSORIE

1.888.229.2875 www.innoviamedical.com/instrusafe Summit Medical
An Innovia Medical Company

711 L x 279 W x 76 H (mm)

Container Foot

IN-8920-TF

InstruSafe Sterile Processing



Preconfigured Trays for da Vinci Xi **Single-Site**



Preconfigured Trays for da Vinci Xi®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Instruments—Single or Dual Staplers and their accommodating accessories—our trays are the recommended solutions for instrument reprocessing.

INSTRUMENT TRAY

Six Instrument Tray	
Holds:	
• 6 (5mm) Xi Single-Site Ins	truments
26 ¾ L x 10 W x 3 ¼ H (in)	679 L x 254 W x 92 H (mm)
Wrapped Foot	Container Foot
IN-8914-R	N/A

All da Vinci Xi Single-Site Trays have green covers.

Excludes IN-8911



☐ INSTRUMENT TRAY


Single Instrument Tray Holds: • 1 EndoWrist Instrument 689 L x 89 W x 54 H (mm) 27 ½ L x 3 ½ W x 2 ½ H (in) Wrapped Foot Container Foot IN-8931-R N/A

INSTRUMENT TRAY

Two Instrument Tray

• 2 EndoWrist instruments

27 1/8 L x 3 3/8 W x 3 1/2 H (in) 689 L x 86 W x 89 H (mm) Container Foot Wrapped Foot IN-8947-R N/A



INSTRUMENT TRAY

Three Instrument Tray

• 3 EndoWrist instruments

26 ³/₄ L x 6 ⁷/₈ W x 3 ³/₄ H (in) 679 L x 175 W x 94 H (mm) Wrapped Foot Container Foot IN-8948-R IN-8948-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing



Preconfigured Trays for

da Vinci Xi



Preconfigured Trays for

da Vinci Xi



INSTRUMENT TRAY

Four Instrument Tray Holds:	
4 EndoWrist Instruments	
26 ¾ L x 6 ½ W x 3 ¾ H (in)	679 L x 175 W x 94 H (mm)
Wrapped Foot IN-8950-R	Container Foot IN-8950-TF



INSTRUMENT TRAY

Three Instrument + Accessory Tray

• 3 EndoWrist Instruments + Associated Accessories

26 ¾ L x 6 ½ W x 3 ¾ H (in)		679 L x 175 W x 94 H (mm)		
	Wrapped Foot IN-8949-R	Container Foot IN-8949-TF		

INSTRUMENT TRAY

Six Instrument Tray Holds: • 6 EndoWrist Instruments	
26 ¾ L x 10 W x 3 ¾ H (in)	679 L x 254 W x 92 H (mm)
Wrapped Foot IN-8936-R	Container Foot N/A



PROCEDURE TRAY

Four Instrument Procedure Tray

• 4 EndoWrist Instruments + Associated Accessories

28 L x 8 W x 5 7/16 H (in)	711 L x 203 W x 138 H (mm)
Wrapped Foot	Container Foot
IN-8915-R	N/A

INSTRUMENT TRAY

Eight Instrument Tray Holds: • 8 EndoWrist Instruments	
26 ¾ L x 13 W x 3 % H (in)	679 L x 330 W x 92 H (mm)
Wrapped Foot IN-8938-R	Container Foot N/A



PROCEDURE TRAY

Six Instrument Procedure Tray Holds:		
26 ¾ L x 13 W x 3 ¾ H (in) 679 L x 330 W x 95 H (mr		
Wrapped Foot	Container Foot	
IN-8937-R	N/A	

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

64 1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing



Preconfigured Trays for

da Vinci Xi



Preconfigured Trays for

da Vinci Xi



PROCEDURE TRAY

Eight Instrument Procedure Tray Holds:			
• 8 EndoWrist Instruments + Associated Accessories			
26 ¾ L x 13 W x 3 % H (in)	679 L x 330 W x 92 H (mm)		
Wrapped Foot IN-8935-R	Container Foot N/A		







ACCESSORY TRAY

Ten	A	CC	es	sor	y T	ray

Holds:

- 6 Cannulas or Obturators
- 1 Release Tool
- 3 Cables

19 % L x 9 % W x 2 H (in)	498 L x 244 W x 51 H (mm)		
Wrapped Foot	Container Foot		
IN-8939-R	IN-8939-TF		

SCOPE TRAY

Single Scope Tray Holds: • 1 (8mm) 0° or 30° Scope	
24 % L x 13 ½ W x 4 ¾ H (in)	632 L x 343 W x 111 H (mm)
Wrapped Foot INS-8930-R	Container Foot N/A

□ CAMERA TRAY

Single Camera Tray	
Holds:	
• 1 Hand Held Camera	
15 ¼ L x 9 5% W x 3 H (in)	387 L x 244 W x 76 H (mm)
Wrapped Foot INS-8934-R	Container Foot INS-8934-TF

NEED A PLACE TO SOAK? **SEE PAGE 90 FOR DETAILS**

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



STAPLER + ACCESSORY TRAY

STAPLER + ACCESSORY TRAY

• 1 Stapler Instrument + Associated Accessories

Single Stapler + Accessory Tray

27 ½ L x 3 ¾ W x 3 ½ H (in)

Wrapped Foot

IN-8933-R

Dual Stapler + Accessory Tray		
Holds:		
• 2 Stapler Instruments + Associated Accessories		
26 ³ / ₄ L x 6 ⁷ / ₈ W x 3 ¹ / ₁₆ H (in)	679 L x 175 W x 94 H (mm)	
Wrapped Foot	Container Foot	
IN-8932-R	N/A	

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

689 L x 86 W x 89 H (mm)

Container Foot

N/A



Preconfigured Trays for da Vinci Si® Single-Site®

Coordinate your instrument sets with multi-level trays that accommodate a full Single-Site instrument and accessory procedure set all the way down to one Single-Site instrument.



Preconfigured Trays for

da Vinci Si Single-Site







□ INSTRUMENT TRAY

Single Instrument Tray Holds:	
• 1 (5mm) Si Single-Site Instrument	
25 % L x 3 % W x 2 % H (in)	651 L x 79 W x 60 H (mm)
Wrapped Foot IN-8901-R	Container Foot N/A

PROCEDURE TRAY

Seven Instrument + Accessory Tray

- 7 (5mm) Si Single-Site Instruments
- Flexible Obturators
- 1 Single-Site Suction Irrigator
- Curved Cannulas

Alignment Target	
22 ½ L x 13 W x 4 ½ H (in)	562 L x 330 H x 105 H (mm)
Wrapped Foot	Container Foot
IN-8902-G2-R	IN-8902-G2-TF
	22 ½ L x 13 W x 4 ½ H (in)

PROCEDURE TRAY

Ten Instrument + Accessory Tray

- 10 (5mm) Si Single-Site Instruments
- Alignment Target
- Flexible Obturators
- Curved Cannulas
- 1 Single-Site Suction Irrigator

23 % L x 9 % W x 6 H (in)	606 L x 244 W x 152 H (mm)
Wrapped Foot	Container Foot
IN-8904-R	IN-8904-TF

All da Vinci Si Single-Site Trays have green covers.

Excludes IN-8901



ACCESSORY TRAY

Accessory Tray

Holds:

- Alignment Target
- Curved Cannulas

Flexible Obturators	
22 ½ L x 13 W x 2 ¾ H (in)	562 L x 330 W x 60 H (mm)
Wrapped Foot IN-8903-R	Container Foot IN-8903-TF

INSTRUMENT + SUCTION TRAY

Seven Instrument + Suction Tray

Holds:

- 7 (5mm) Si Single-Site Instruments
- 1 Single-Site Suction Irrigator

22 ½ L x 13 W x 2 3/8 H (in)	512 L x 330 W x 60 H (mm)
Wrapped Foot IN-8907-R	Container Foot IN-8907-TF

All da Vinci Si Single-Site Trays have green covers.

Excludes IN-8901

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing



Preconfigured Trays for da Vinci Si **EndoWrist**®

No two sets are alike. That is why we have created a wide range of trays to accommodate single or multiple EndoWrist instruments and their accompanying accessories for the Si Surgical System.




Preconfigured Trays for

da Vinci Si **EndoWrist**









INSTRUMENT TRAY

Single Instrument Tray Holds: • 1 EndoWrist Instrument	
25 ¼ L x 3 ¼ W x 2 ¾ H (in)	641 L x 83 W x 60 H (mm)
Wrapped Foot IN-8885-R	Container Foot N/A



□ INSTRUMENT TRAY

Ten Instrument Tray	
Holds: • 10 EndoWrist Instruments	
23 % L x 9 % W x 4 % H (in)	606 L x 244 W x 105 H (mm)
Wrapped Foot IN-8898-01-R	Container Foot IN-8898-01-TF



□ INSTRUMENT TRAY

Five Instrument Tray Holds: • 5 EndoWrist Instruments	
23 % L x 9 % W x 2 % H (in)	606 L x 244 W x 60 H (mm)
Wrapped Foot IN-8894-R	Container Foot IN-8894-TF



☐ INSTRUMENT TRAY

Twelve Instrument Tray Holds: • 12 EndoWrist Instruments	
22 ½ L x 13 W x 4 ½ H (in)	572 L x 330 W x 114 H (mm)
Wrapped Foot IN-8892-R	Container Foot IN-8892-TF



INSTRUMENT TRAY

Six Instrument Tray Holds: • 6 EndoWrist Instruments	
22 ½ L x 13 W x 2 ¾ H (in)	562 L x 330 W x 60 H (mm)
Wrapped Foot IN-9999-160-01-R	Container Foot IN-9999-160-01-TF



PROCEDURE TRAY

PROCEDURE TRAY	
Eight Instrument Proced	lure Tray
Holds:	
8 EndoWrist Instruments	
• 1 Alignment Target	
• 3 Obturators	
• 3 Cannulas	
23 % L x 9 % W x 4 % H (in)	606 L x 244 W x 105 H (mm)
Wrapped Foot	Container Foot
IN-8984-01-R	IN-8984-01-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

1.888.229.2875 70 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing

ROBOTICS ROBOTICS



Preconfigured Trays for

da Vinci Si **EndoWrist**



Preconfigured Trays for

da Vinci Si® Stapler

Get more flexibility to meet the demands of your procedure schedule with trays designed specifically for use with the Si Series da Vinci Surgical Stapler System. Options include various feet choices for use with wrap or rigid containers.



PROCEDURE TRAY

Ten Instrument Procedure Tray

Holds:

- 10 EndoWrist Instruments
- 1 Alignment Target
- 3 Obturators
- 3 Cannulas

3 Callilulas	
22 ½ L x 13 W x 4 ½ H (in)	572 L x 330 W x 114 H (mm)
Wrapped Foot IN-8982-R	Container Foot IN-8982-TF



TAPLER TRAY


Holds:

- 1 EndoWrist Stapler Instrument
- 1 EndoWrist 13mm Stapler Cannula

Single Stapler + Accessory Tray

- 1 EndoWrist 13mm Stapler Obturator • 1 EndoWrist 13-8mm Stapler Reducer
- 1 EndoWrist Stapler Cable
- 1 EndoWrist Stapler Release Kit
- 1 Instruction Plaque

25 L x 5 3/8 W x 3 H (in)	635 L x 137 W x 76 H (mm)
Wrapped Foot	Container Foot
IN-8940-R	IN-8940-TF



ACCESSORY TRAY

Eight Accessory Tray

- 1 Alignment Target
- 3 Obturators
- 1-2 Cables
- 3 Cannulas

o daminatas	
15 ¼ L x 9 5/8 W x 3 H (in)	387 L x 244 W x 76 H (mm)
Wrapped Foot IN-8899-R	Container Foot IN-8899-TF



STAPLER TRAY

Single Stapler + Accessory Tray

Holds:

- 1 EndoWrist Stapler Instrument
- 1 EndoWrist 13mm Stapler Cannula
- 1 EndoWrist 13mm Stapler Obturator
- 1 EndoWrist 13-8mm Stapler Reducer
- 1 EndoWrist Stapler Cable
- 1 EndoWrist Stapler Release Kit
- 1 Instruction Plaque

23 % L x 9 % W x 2 % H (in) 606 Wrapped Foot	6 L x 244 W x 60 H (mm)
Wrapped Foot	
IN-8944-R	Container Foot IN-8944-TF




ACCESSORY TRAY

Seven Accessory Tray

Holds:

- 1 Alignment Target
- 3 Obturators
- 3 Cannulas

• 5 Cannulas	
8 % L x 9 % W x 3 H (in)	213 L x 244 W x 76 H (mm)
Wrapped Foot IN-8897-R	Container Foot IN-8897-TF



₽ STAPLER TRAY

Dual Stapler Instrument Tray	
Holds:	
• 1-2 EndoWrist Stapler Instruments	
25 L x 5 % W x 2 % H (in)	635 L x 137 W x 60 H (mm)
Wrapped Foot	Container Foot
IN-8945-R	N/A

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing

ROBOTICS ROBOTICS



Preconfigured Trays for

da Vinci Si Stapler




Preconfigured Trays for da Vinci Si Scope

One of our original da Vinci products, Single and Dual Scope Trays are designed to protect the full variety of available endoscopes for the Si Surgical System and are 510(k) cleared for use in the STERRAD® 100NX® Express cycle.











STAPLER TRAY

Dual Stapler + Accessory Tray

- 2 EndoWrist Stapler Instruments
- 2 EndoWrist 13mm Stapler Cannulas
- 2 EndoWrist 13mm Stapler Obturators
- 2 EndoWrist 13-8mm Stapler Reducers
- 2 EndoWrist Stapler Cables
- 1 EndoWrist Stapler Release Kit
- 1 Instruction Plaque

23 % L x 9 % W x 2 % H (in)	606 L x 244 W x 60 H (mm)
Wrapped Foot	Container Foot
IN-8942-R	IN-8942-TF



Single Si Stapler Motor Pack Accessory Tray Holds: • 1 EndoWrist Stapler 45 Motor Pack			
		8 ¼ L x 9 % W x 3 H (in)	210 L x 244 W x 76 H (mm)
		Wrapped Foot IN-8943-R	Container Foot IN-8943-TF



SCOPE TRAY

Single Scope Tray • 1 (12mm) 0° or 30° Scope 22 L x 5 5/8 W x 3 H (in) 559 L x 143 W x 76 H (mm) Wrapped Foot Container Foot IN-8891-SI-12-S-R N/A

SCOPE TRAY

Single Scope Tray	
Holds:	
• 1 (8.5mm) 0° or 30° Scope	
22 L x 5 5 W x 3 H (in)	559 L x 143 W x 76 H (mm)
Wrapped Foot IN-8891–SI-85–S-R	Container Foot N/A

ACCESSORY TRAY

Stapler Accessory Tray

- 2 EndoWrist Stapler Cannulas
- 2 EndoWrist Stapler Cannula Reducers
- 2 EndoWrist Stapler Cannula Blunt Obturators • 2 EndoWrist Stapler Cables
- 1 EndoWrist Stapler Release Kit

25 L x 5 3/8 W x 2 3/8 H (in)	635 L x 137 W x 60 H (mm)
Wrapped Foot	Container Foot
IN-8946-R	IN-8946-TF



SCOPE TRAY

Dual Scope Tray	
Holds:	
• 2 (12mm) 0° or 30° Scope:	S
25 L x 5 5 W x 3 H (in)	635 L x 143 W x 76 H (mm)
Wrapped Foot IN-9999-172-S-R	Container Foot N/A

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing

ROBOTICS ROBOTICS



Preconfigured Trays for

da Vinci Si Scope



Preconfigured Trays for

Senhance[™] Surgical **Robotic Trays**

A sophisticated robotic system deserves a solution just as sophisticated to hold its instruments. That's why TransEnterix collaborated with us, the instrument protection tray experts. This line of trays is not only designed and constructed to fit the Senhance Surgical Robotic System instruments and adapters like a glove, but it is also validated specifically for this instrumentation.



SCOPE TRAY

Dual Scope Tray Holds: • 2 (8.5mm) 0° or 30° Scope	e
25 L x 5 % W x 3 H (in)	635 L x 143 W x 76 H (mm)
Wrapped Foot IN-9999-178-S-R	Container Foot N/A




INSTRUMENT TRAY

Instrument Tray

Holds:

- 8 Instrument Adapters
- 12 Instrument Shafts
- Cables

19 5% L x 9 5% W x 4 1% H (in)	498 L x 244 W x 105 H (mm)
Wrapped Foot	Container Foot
IN-8750-R	IN-8750-TF



SCOPE TRAY

Dual Scope Tray	
Holds:	
• 2 (12mm) 0° or 30° Scopes	
21 ¼ L x 9 W x 3 H (in)	540 L x 229 W x 76 H (mm)
Wrapped Foot IN-8989-S-R	Container Foot IN-8989-S-TF



ULTRASONIC INSTRUMENT TRAY

Instrument Tray

Holds:

- 1 Ultrasonic Adapter
- 1 Ultrasonic Transducer

TOTTI ASOTTIC TTATISAUCCI	
25 L x 3 3/8 W x 3 H (in)	635 L x 137 W x 76 H (mm)
Wrapped Foot IN-8753-R	Container Foot IN-8753-TF









SCOPE TRAY

Dual Scope Tray	
Holds:	
• 2 (8.5mm) 0° or 30° Scop	es
21 ¼ L x 9 W x 3 H (in)	540 L x 229 W x 76 H (mm)
Wrapped Foot	Container Foot
IN-8988-S-R	IN-8988-S-TF





ENDOSCOPE ADAPTER TRAY

Single Instrument Tray	
Holds:	
• 1 Endoscope Adapter	
19 % L x 9 % W x 2 % H (in)	498 L x 244 W x 60 H (mm)
Wrapped Foot IN-8751-R	Container Foot IN-8751-TF

ENDOSCOPE ADAPTER TRAY

Dual Instrument Tray Holds:

• 2 Endoscope Adapters


l	2 Elidoscope Adapters	
	19 ½ L x 9 ½ W x 2 ½ H (in)	498 L x 244 W x 60 H (mm)
	Wrapped Foot IN-8752-R	Container Foot IN-8752-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

Senhance is a trademark of TransEnterix Surgical, Inc. Please refer to the Senhance™ Surgical Robotic System IFU for instrument reprocessing.

1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing



Custom Configured Tray Example for

Robotic Instrumentation



Preconfigured Trays for Rigid Scopes

Our rigid scope trays are universal and designed to fit

ROBOTIC LAPAROSCOPIC SET | 9 INSTRUMENTS

CUSTOM CONFIGURED

AFTER

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc. Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection



scopes of any size for any procedure.



RIGID SCOPE TRAY 4-7"

• Up to 3 (2-10mm) Scopes

• Accommodates scopes from 4-7" in working length

11 L x 3 ¼ W x 2 ¾ H (in) 279 L x 83 W x 60 H (mm) IN-6011-R N/A



RIGID SCOPE TRAY 5-131/2"

Holds:

• Up to 3 (2-10mm) Scopes

• Accommodates scopes from 5-13½" in working length

17 ½ L x 3 ¼ W x 2 ¾ H (in) 445 L x 83 W x 60 H (mm) Container Foot IN-6012-R N/A



RIGID SCOPE TRAY 7-161/2"

• Up to 3 (2-10mm) Scopes

• Accommodates scopes from 7-16½" in working length

20 ½ L x 3 ½ W x 2 ¾ H (in) 515 L x 83 W x 60 H (mm) Wrapped Foot IN-6013-R N/A

SCOPES SCOPES



Preconfigured Trays for

Rigid Scopes



Preconfigured Trays for Scope + Cable

Our scope + cable trays allow for faster loading before sterilization while also protecting the delicate fiber optics within.

RIGID SCOPE TRAY 91/2-211/2"

• Up to 3 (2-10mm) Scopes

• Accommodates scopes from 9½-21½" in working

25 ½ L x 3 ½ W x 2 ¾ H (in)	641 L x 83 W x 60 H (mm)
Wrapped Foot IN-6014-R	Container Foot N/A



• Up to 3 (2-10mm) Scopes

• Accommodates scopes from 11-26" in working length

30 ½ L x 3 ¼ W x 2 ¾ H (in)	765 L x 83 W x 60 H (mm)
Wrapped Foot	Container Foot
IN-6015-R	N/A



SCOPE + CABLE TRAY 5-6"

Holds:

• Up to 3 (2-10mm) Scopes

• 1 Cable

• Accommodates scopes from 5-6" in working length

1	11 L x 6 3/8 W x 2 3/8 H (in)	279 L x 162 W x 60 H (mm)
	Wrapped Foot	Container Foot
	IN-6031-R	N/A

SCOPE + CABLE TRAY 5-131/2"

Holds:

• Up to 3 (2-10mm) Scopes

• 1 Cable

• Accommodates scopes from 5-13½" in working length

20 ½ L x 5 5% W x 2 3% H (in)	521 L x 143 W x 60 H (mm)
Wrapped Foot IN-6032-R	Container Foot N/A

SCOPE + CABLE TRAY 6½-17"

Holds:


• Up to 3 (2-10mm) Scopes

• 1 Cable

• Accommodates scopes from 6½-17" in working length

25 L x 5 3/8 W x 2 3/8 H (in)	635 L x 137 W x 60 H (mm)
Wrapped Foot	Container Foot
IN-6033-R	N/A

SCOPES



SCOPES

Custom Configured Tray Example for

Scope + Cable

BEFORE



Preconfigured Trays for Flexible Scopes

Our durable flexible scope trays are sure to protect your delicate and expensive scopes. Our Small Parts Boxes are perfect to keep all your accessories

RESECTOSCOPE | 11 INSTRUMENTS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





FLEXIBLE SCOPE TRAY

Holds: • 1 Flexible Scope • Accommodates scopes from 12-27" in working length 19 % L x 9 % W x 3 H (in) 498 L x 244 W x 76 H (mm) Wrapped Foot IN-6021-R N/A

FLEXIBLE SCOPE TRAY

Holds:	
• 1 Flexible Scope with Atta	ched Light Source
• Accommodates scopes fro	om 12-27" in working length
19 % L x 9 % W x 3 H (in)	498 L x 244 W x 76 H (mm)
Wrapped Foot	Container Foot
IN-6022-R	N/A

FLEXIBLE SCOPE TRAY

• 1 Flexible Video Scope	
• Accommodates scopes fro	om 12-27" in working length
19 % L x 9 % W x 3 H (in)	498 L x 244 W x 76 H (mm
Wrapped Foot	Container Foot
IN-6023-R	N/A


1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing



Preconfigured Tray for

Flexible Scopes



Custom Configured Tray Example for

Flexible Scopes




FLEXIBLE SCOPE TRAY


Holds:		
• 1 Flexible Bronchscope		
• Accommodates scopes from 12-27" in working length		
19 % L x 9 % W x 3 H (in)	498 L x 244 W x 76 H (mm)	
Wrapped Foot IN-6024-R	Container Foot N/A	

FLEXIBLE URETEROSCOPE | 2 INSTRUMENTS

BEFORE




 $\textit{Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection and the provide 360 degrees of provide 360 degrees of provide 360 degrees of protection and the provide 360 degrees of provide 360 degrees of protection and the protection and the provide 360 degrees of protection and the protection and the protection and the provide 360 degrees of protection and the pro$



Custom Configured Tray Example for

Flexible Scopes



Custom Configured Tray Example for

Flexible Scopes


VIDEO LARYNGOSCOPE | 1 INSTRUMENT

BEFORE

COLONOSCOPY TRAY | 1 INSTRUMENT




AFTER



Custom Configured Tray Example for

Flexible Scopes




Custom Configured Tray Example for

Flexible Scopes

FLEXIBLE SCOPE TRAY | 3 INSTRUMENTS

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

FLEXIBLE SCOPE TRAY | 1 INSTRUMENT



88 1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing

TRANSPORT CONTAINERS TRANSPORT CONTAINERS

XL Transport Containers

Transporting da Vinci SP, Xi scopes, and EndoWrist instruments appropriately and effectively is a difficult task. Use the InstruSafe XL Transport Container to get the job done right! This durable container is the perfect size, enabling you to soak and protect your instruments while they travel from the O.R. to sterile processing. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.

XL Transport Container

Special Features

- Highly durable material designed to withstand everyday use
- Cart wash safe—maximum temperature of 230°F
- Economical and creates a dedicated space for horizontally soaking robotic instruments
- Made in the USA



XL TRANSPORT CONTAINER

XL Transport Container + Optional Drain

• 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments + Associated Accessories

34 L x 17 ¼ W x 6 H (in) 864 L x 438 W x 152 H (mm)

With Drain

Without Drain

IN-PIC-20-D

IN-PIC-20

IN-PIC-20-D-LT

XL TRANSPORT CONTAINER SYSTEM

XL Transport Container with Lifting Tray + Optional Drain

- Lifting Tray
- 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments + Associated Accessories

34 L x 17 1/4 W x 6 H (in) 864 L x 438 W x 152 H (mm)

With Drain IN-PIC-20-D-LT IN-PIC-20-LT



LIFTING TRAY

Single Endoscope Configuration

• 1 da Vinci Xi or SP Endoscope + Associated Accessories Multiple Instrument Configuration

• Up to 8 da Vinci Xi or SP EndoWrist Instruments + Associated Accessories

27 ½ L x 13 ½ W x 3 % H (in)

699 L x 343 W x 98 H (mm)

IN-PIC-34





Holds: up to 8 da Vinci Xi or SP EndoWrist Instruments



Heavy Duty Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. They are similar to our other transport containers but have a more robust construction and feature stainless steel hinges and latches that can be easily undone with just one hand.


Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. Also, you can easily identify content through the transparent amber color. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.



HEAVY DUTY TRANSPORT CONTAINER

Short Container	
20 ½ L x 12 5/8 W x 3 3/4 H (in)	521 L x 321 W x 95 H (mm)
IN-PIC-15	



Tall Container	
20 L x 12 ¼ W x 6 ½ H (in)	510 L x 310 W x 156 H (mm)
IN-PIC-16	



Tall Container	
20 L x 12 ¼ W x 6 ½ H (in)	510 L x 310 W x 156 H (mm)
IN-PIC-16	

SPECIAL FEATURES

- Autoclavable
- Made of PEI material (high-temperature plastic)
- High durability against automated washing
- Transparent amber color for easy instrument identification
- Easy locking, removable cover with removable silicone gasket that can be replaced
- Lightweight
- Combine with Instru-Zyme for optimal pre-cleaning



TRANSPORT CONTAINER

	Small, Short Container		
13 % 16 L x 15 5/8 W x 4 5/16 H (in)		345 L x 398 W x 110 H (mm)	
	IN-PIC-14		

TRANSPORT CONTAINER

Large, Short Container	
20 % L x 15 % W x 5 5/16 H (in)	524 L x 397 W x 135 H (mm)
IN-PIC-12	

TRANSPORT CONTAINER

Small, Tall Container	
13 % ₁₆ L x 15 % W x 7 ½ H (in)	345 L x 398 W x 190 H (mm)
IN-PIC-13	

TRANSPORT CONTAINER

Large, Tall Container		
20 % L x 15 % W x 7 ½ H (in)	524 L x 397 W x 191 H (mm)	
IN-PIC-11		

1.888.229.2875 www.innoviamedical.com/instrusafe

Summit Medical An Innovia Medical Company InstruSafe Sterile Processing

ACCESSORIES ACCESSORIES

Silicone Small Parts Boxes

Don't sweat the small stuff. Made of silicone, these boxes are perfect for keeping extra delicate small parts and accessories organized and protected.

Plaque Options


A variety of plaque colors are available to help identify specific sets by procedure or service type. Affix them to our trays and modules using our standard screws. Common identification uses are for bar code placement, tray name labels, or other custom information. Laser engraving of plaques is provided at an additional cost.

ARIATRIC

SILICONE SMALL PARTS BOX

3x3" Small Parts Box	
3 L x 3 W x ¾ H (in)	76 L x 76 W x 19 H (mm)
IN-2652	

Our in-house laser engraving capabilities allow us to personalize trays with custom identification and/ or specifications, such as interior instrument descriptions, illustrations, exterior procedure types, specialty, or facility information.



SILICONE SMALL PARTS BOX

3x6" Small Parts Box		
6 L x 3 W x ³ / ₄ H (in)	152 L x 76 W x 19 H (mm)	
IN-2653		

PLAQUES

2 L x 2 1/8 W (in)

51 L x 56 W (mm)

Plaque Assembly for Trays		
1 L x 4 ¾ W (in) 25 L x 111 (mm)	IN-2550-01	
1 L x 4 3/8 W (in) 25 L x 111 (mm)	IN-2550-05	BLACK
1 L x 4 ¾ W (in) 25 L x 111 (mm)	IN-2550-15	BLACK
Please indicate desired color of plaque using corresponding number below in place of XX using the color options table.		
4 3/8 L x 2 W (in) 111 L x 51 W (mm)	IN-2600-XX •••	

IN-2300-XX •••

000000

Steam Sterilization
These anodized aluminum plaques for use with

**COLOR OPTIONS

**COLOR OPTIONS

**COLOR OPTIONS

**Low-Temp Sterilization

Plaques have a coating that specially withstands low-specially with low-specially with low-specially with low-specially with low

Plaques have a coating that specially withstands low-temperature sterilization. steam only. 01 SILVER N/A 02 RED 12 RED 03 GREEN GREEN 04 14 BLUE BLUE 05 BLACK **BLACK** 06 GOLD GOLD 07 PURPLE PURPLE 08 MINT GREEN

ROBOTICS

ANSPORT CONTAINERS


Made entirely of medical grade silicone

• Maximum protection for small parts

• Tethered cover to prevent accidental loss

Can be attached to trays or placed on posts for easy removal

• Textured internal surface and proper ventilation for efficient sterilant penetration



Day Ear Hook 2mm

Buck Ear Curette #1
straight

Buck Ear Curette #2
Straight

Buck Ear Curette #1
Curved

Buck Ear Curette #2
Curved

Ear Knife Sickle Fine

Ear Knife 1mm Round

Billeau Ear Loop 2.5mm
Small

Shapleigh Ear Curette
3.5mm

ACCESSORIES ACCESSORIES

Finger Mats

These textured silicone instrument cushion mats are available in a variety of sizes and textures to fit all of your instrument protection needs.



FINE FINGER MAT

Full Sheet		
19 L x 13 W (in)	483 L x 330 W (mm)	
IN-2616		

FINE FINGER MAT

Half Sheet	
9 ½ L x 13 W (in)	241 L x 330 W (mm)
IN-2617	

FINE FINGER MAT

Quarter Sheet	
4 ³ / ₄ L x 13 W (in)	121 L x 330 W (mm)
IN-2	2618

Guardian™ Tray

Simple, small instrument sets now have a new home. Our silicone finger mats with sidewalls provide a safe and pliable environment for them to reside in.



GUARDIAN TRAY

Anodized Aluminum Tray with Finger Mat		
10 L x 6 W x 1 ¼ H (in)	254 L x 152 W x 32 H (mm)	
IN-1310		

- · Accommodate a wide variety of small instrumentation
- Feature built-in and vented silicone finger mats
- Reduce instrument repair and replacement costs
- Allow maximum sterilant penetration
- Prevent moisture accumulation

Instru-Zyme[®]

This powerful enzyme and detergent formulation is designed to assist in the removal of blood, tissue, mucus and other proteinaceous materials. Use this product with our transport containers for a simple solution to water and detergent.



INSTRU-ZYME

• Specific Gravity: 1.200 • Color: Light Blue • pH: 7.8

• Tropical Fragrance

• Form: Liquid Gel

IN-1032

• Will not stain instruments or equipment made of: metal, plastic, glass or rubber

• Non-corrosive

Biodegradable

6 Bottles/Case

• Sustains moisture during transport

• Provides a ready-to-use, neutral pH enzyme solution

• Reduces heavy instrument transport weight as compared to liquid-filled trays

• Minimizes potential exposure due to splashing

• Spray applicator for even coverage of soiled areas

No rinsing necessary before automated processing

Rounders® Point Guards

Wrap your trays with confidence. Point Guards are engineered to sit securely on almost any tray corner. Prevent holes and tears in your sterilization wrap caused by surgical trays that are heavy, have sharp corners or are improperly handled.



ROUNDERS POINT GUARDS

IN-7330 80/Bag

Air cushion tray corner protection

• Help avoid compromised sterility in the O.R.

• Can be used in steam, ethylene oxide and low temperature plasma sterilization

• Reusable and economical

• Patent #9119611

Non-sterile

MADE IN THE USA

WISH LIST

InstruSafe Sterile Processing

Summit Medical LLC An Innovia Medical Company 815 Vikings Parkway Suite 100 St. Paul, MN 55121 | USA

1.888.229.2875 www.innoviamedical.com/instrusafe

InstruSafe Instrument Protection Trays are proudly made in St. Paul, Minnesota, USA.

 $Summit\ Medical\ LLC\ is\ an\ FDA\ registered\ facility, ISO\ 13485: 2016\ certified\ and\ CMDCAS\ certified.$

 $In struSafe\ Trays\ are\ 510(k)\ cleared\ for\ multiple\ sterilization\ cycles\ to\ give\ our\ customers\ more\ options.$

CAUTION: Federal (USA) law restricts these devices for sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

Specifications and designs are subject to change without any notice or obligation on the part of the manufacturer.

da Vinci, SP, S, Si, Xi, Single-Site and EndoWrist are registered trademarks of Intuitive Surgical, Inc. InstruSafe, Instru-Zyme, Hold-Its, Hold-Downs, Rounders Point Guards, and Guardian are trademarks and registered trademarks of Summit Medical.

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.

All brands are trademarks or registered trademarks of their respective owners.

©2020 Innovia Medical. All rights reserved. M-5 | RevB | DEC2020